

Product datasheet

Anti-MCM9 antibody ab103420

2 图像

概述

产品名称	Anti-MCM9抗体
描述	兔多克隆抗体to MCM9
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Rat, Human
免疫原	Full length Human MCM9 Isoform 2 (NP_694987.1, amino acids 1-391).
阳性对照	Rat brain tissue lysate; MCM9 transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab103420** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: Isoform 1, 127 ; Isoform 2, 44 kDa.

靶标

相关性 MCM9, a member of the MCM2-8 family, binds to chromatin and is required for the recruitment


of the MCM2-7 helicase onto chromatin. MCM9 can form a complex with Cdt1. It is thought that MCM9 might play an important role in DNA replication since its depletion results in its inhibition.

图片



Anti-MCM9 antibody (ab103420) at 5 µg/ml +
Rat brain tissue lysate at 50 µg

Predicted band size: Isoform 1, 127 ; Isoform
2, 44 kDa



All lanes : Anti-MCM9 antibody (ab103420)
at 5 µg/ml

Lane 1 : MCM9 transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Predicted band size: Isoform 1, 127 ; Isoform
2, 44 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors