

## Product datasheet

# Anti-MCSF Receptor antibody ab10676

### 概述

产品名称	Anti-MCSF Receptor抗体
描述	山羊多克隆抗体to MCSF Receptor
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Human
免疫原	Recombinant fragment (extracellular domain) (Human).

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab10676** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		
WB		

应用说明	<p>ELISA: Use at a concentration of 0.5 - 1.0 µg/ml. The detection limit for recombinant human M-CSF Receptor is approximately 0.15 ng/well.</p> <p>WB: Use at a concentration of 0.1 - 0.2 µg/ml. The detection limit for recombinant human M-CSF Receptor is approximately 5 ng/lane under non-reducing and 25 ng/lane under reducing conditions. Predicted molecular weight: 81 kDa.</p>
------	---

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

---

## 靶标

<b>功能</b>	Protein tyrosine-kinase transmembrane receptor for CSF1 and IL34.
<b>组织特异性</b>	Expressed in bone marrow and in differentiated blood mononuclear cells.
<b>序列相似性</b>	Belongs to the protein kinase superfamily. Tyr protein kinase family. CSF-1/PDGF receptor subfamily. Contains 5 Ig-like C2-type (immunoglobulin-like) domains. Contains 1 protein kinase domain.
<b>细胞定位</b>	Membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

---

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

---

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors