

Product datasheet

Anti-MGC4618 antibody ab49189

1 图像

概述

产品名称	Anti-MGC4618抗体
描述	兔多克隆抗体to MGC4618
宿主	Rabbit
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Horse, Guinea pig, Cat, Dog
免疫原	A region within synthetic peptide: MLSFSDALLS IIATVMILPV THTESISPEQQ FDRSVQRLLA TRIAVYLMTF, corresponding to N terminal amino acids 36-85 of Human MGC4618 Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	HepG2 cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab49189** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明

应用	Ab 评论	说明
----	-------	----

WB Use a concentration of 1.25 µg/ml. Detects a band of approximately 56 kDa (predicted molecular weight: 56 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.


ELISA Use at an assay dependent concentration.
ELISA titre using peptide based assay: 1:1562500.

靶标

序列相似性 Belongs to the TMEM175 family.

细胞定位 Membrane.

图片



Anti-MGC4618 antibody (ab49189) at 1.25 µg/ml + HepG2 cell lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

Developed using the ECL technique.

Predicted band size: 56 kDa

Observed band size: 56 kDa

Diluted in 5% skim milk / PBS buffer, gel concentration 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors