

Product datasheet

Anti-MIS12 antibody ab92402

1 图像

概述

产品名称	Anti-MIS12抗体
描述	兔多克隆抗体to MIS12
宿主	Rabbit
经测试应用	适用于: IHC-P
种属反应性	与反应: Human 预测可用于: Chimpanzee, Gorilla, Orangutan, Bat
免疫原	Synthetic peptide from within residues: AKLKQTLTFF DELHNVGRDH GTSDFRESLV SLVQNSRKLQ NIRDNVEKES KRLKIS, corresponding to C terminal amino acids 150-205 of Human MIS12 (NP_076944.1). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Human stomach adenocarcinoma tissue

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab92402** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

IHC-P

应用说明

IHC-P: 1/100 - 1/500. Antigen retrieval recommended.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标


相关性

MIS12 is part of a complex that plays an essential role in chromosome segregation in vertebrates and contributes to mitotic kinetochore assembly.

细胞定位

Nucleus. Associated with the kinetochore.

图片



Formalin-fixed, paraffin-embedded section of Human stomach adenocarcinoma stained with ab92402 used at a dilution of 1/250.

Detection: DAB staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-MIS12 antibody (ab92402)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors