

Product datasheet

Anti-Methicillin Resistant Staphylococcus Aureus antibody [MRSA 198] ab73263

1 References

概述

| | |
|-------|---|
| 产品名称 | Anti-Methicillin Resistant Staphylococcus Aureus 抗体 [MRSA 198] |
| 描述 | 小鼠单克隆抗体 [MRSA 198] to Methicillin Resistant Staphylococcus Aureus |
| 特异性 | ab73263 recognizes the product of the MRSA gene expressed by methicillin-resistant Staphylococcus aureus. |
| 经测试应用 | 适用于: WB, ELISA |
| 种属反应性 | 与反应: Staphylococcus aureus |
| 免疫原 | Gel-purified penicillin-binding protein of MRSA. |
| 常规说明 | Abcam is committed to meeting high standards of ethical manufacturing and has decided to discontinue this product by June 2019 as it has been generated by the ascites method. We are sorry for any inconvenience this may cause. |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| 存储溶液 | Preservative: None Constituents: PBS |
| 克隆 | 单克隆 |
| 克隆编号 | MRSA 198 |
| 同种型 | IgG1 |

应用

Our [Abpromise guarantee](#) covers the use of **ab73263** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----|------|----|
| WB | | |

| 应用 | Ab评论 | 说明 |
|----|------|----|
|----|------|----|

ELISA

| | |
|------|---|
| 应用说明 | <p>ELISA: 1/5000 against purified MRSA protein or MRSA-positive Staphylococcus aureus lysate. WB: 1/2000. Detects a band of approximately 70 kDa.</p> <p>Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the end user.</p> |
|------|---|

靶标

| | |
|-----|---|
| 相关性 | <p>Methicillin-resistant Staphylococcus Aureus (MRSA) is a type of bacteria that is resistant to certain antibiotics. These antibiotics include methicillin and other more common antibiotics such as oxacillin, penicillin and amoxicillin. Staphylococcus infections, including MRSA, occur most frequently among persons in hospitals and healthcare facilities (such as nursing homes and dialysis centres) who have weakened immune systems.</p> |
|-----|---|

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors