

## Product datasheet

# Mouse monoclonal IgG1 [B11/6] (OC515) - Isotype Control ab126030

### 概述

产品名称	小鼠单克隆抗体IgG1 [B11/6] (OC515) -同型对照
偶联物	OC515. Ex: 405nm, Em: 515nm
共轭说明	Abs/Em Max: 405/518 nm
经测试应用	适用于: Flow Cyt
免疫原	Synthetic hapten, which is normally not present in humans or animals.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	pH: 7.40 Preservative: 0.09% Sodium azide Constituent: PBS Note: Contains Antibody Stabiliser.
纯度	Immunogen affinity purified
克隆	单克隆
克隆编号	B11/6
同种型	IgG1
轻链类型	kappa

### 应用

Our [Abpromise guarantee](#) covers the use of **ab126030** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Flow Cyt		Use 5µl for 10 <sup>6</sup> cells.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors