

Product datasheet

Mouse monoclonal [MARG2a-1] Anti-Rat IgG2a heavy chain (FITC) ab11671

1 References

概述

产品名称	小鼠单克隆抗体[MARG2a-1]抗大鼠IgG2a heavy chain (FITC)
宿主	Mouse
靶标种属	Rat
特异性	This antibody specifically recognises rat IgG2a, no cross reactivity is observed with other rat immunoglobulin subclasses.
经测试应用	适用于: Flow Cyt
免疫原	Purified rat immunoglobulins from the IR33 - anti horse ferritin antibody.
偶联物	FITC. Ex: 493nm, Em: 528nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: 50% Glycerol, PBS
纯度	Protein G purified
纯化说明	Purification by affinity chromatography.
克隆	单克隆
克隆编号	MARG2a-1
同种型	IgG1
常规说明	This product should be stored undiluted. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

应用

Our [Abpromise guarantee](#) covers the use of **ab11671** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Flow Cyt		
应用说明	<p>Flow Cyt: Use at a concentration of 5 µg/ml. Use 10µl of the suggested working dilution to label 10⁶ cells in 100µl.</p> <p>Not tested in other applications. Optimal dilutions/concentrations should be determined by the end user.</p>	
<p>Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"</p>		

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors