

Product datasheet

Anti-NABC1 antibody ab106661

3 图像

概述

产品名称	Anti-NABC1抗体
描述	兔多克隆抗体to NABC1
宿主	Rabbit
经测试应用	适用于: IHC-P, ICC/IF, WB
种属反应性	与反应: Human
免疫原	A 16 amino acid synthetic peptide near the C terminal of Human NABC1 (EAW75575).
阳性对照	Human lung tissue lysate; Human breast carcinoma tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (stable for up to 12 months).
存储溶液	Preservative: 0.02% Sodium Azide Constituents: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab106661** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		Use a concentration of 5 µg/ml.
ICC/IF		Use a concentration of 20 µg/ml.
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 62 kDa.

靶标


组织特异性

Expressed in brain and prostate, and at lower levels in testis, intestine and colon.
Overexpressed in most breast cancer cell lines and down-regulated in some colorectal tumors.

细胞定位

Cytoplasm.


图片



Anti-NABC1 antibody (ab106661) at 1 μ g/ml +
Human lung tissue lysate at 15 μ g


Predicted band size: 62 kDa

Western blot - Anti-NABC1 antibody (ab106661)



ab106661, at 5 μ g/ml, staining NABC1 in
Human breast carcinoma tissue by
Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-
embedded sections) - Anti-NABC1 antibody
(ab106661)



Immunofluorescence of NABC1 in Human
Breast Carcinoma cells using ab106661 at 20
 μ g/ml.

Immunocytochemistry/ Immunofluorescence - Anti-
NABC1 antibody (ab106661)

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors