

## Product datasheet

# Anti-NALP12 antibody - Carboxyterminal end ab64928

1 References 1 图像

### 概述

产品名称	Anti-NALP12抗体- Carboxyterminal end
描述	兔多克隆抗体to NALP12 - Carboxyterminal end
经测试应用	适用于: ELISA, IHC-P
种属反应性	与反应: Human
免疫原	This information is considered to be commercially sensitive
阳性对照	Human prostate carcinoma tissue.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Ammonium Sulphate Precipitation
纯化说明	ab64928 is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab64928** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		1/1000.
IHC-P		1/10 - 1/50.

## 靶标

<b>功能</b>	May mediate activation of CASP1 via ASC and promote activation of NF-kappa-B via IKK.
<b>组织特异性</b>	Detected only in peripheral blood leukocytes, predominantly in eosinophils and granulocytes, and at lower levels in monocytes.
<b>疾病相关</b>	Defects in NLRP12 are the cause of familial cold autoinflammatory syndrome type 2 (FCAS2) [MIM:611762]. FCAS are rare autosomal dominant systemic inflammatory diseases characterized by episodes of rash, arthralgia, fever and conjunctivitis after generalized exposure to cold.
<b>序列相似性</b>	Belongs to the NLRP family. Contains 1 DAPIN domain. Contains 8 LRR (leucine-rich) repeats. Contains 1 NACHT domain.
<b>细胞定位</b>	Cytoplasm.

## 图片



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - NALP12 antibody (ab64928)

ab64928 at 1/10 dilution, staining NALP12 in human prostate carcinoma by Immunohistochemistry, Formalin-fixed, Paraffin-embedded tissue, followed by peroxidase-conjugated secondary antibody and DAB staining.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors