

Product datasheet

Anti-NCKAP1L antibody ab82664

1 Abreviews 1 图像

概述

| | |
|-------|---|
| 产品名称 | Anti-NCKAP1L抗体 |
| 描述 | 兔多克隆抗体to NCKAP1L |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: WB, ELISA |
| 种属反应性 | 与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Guinea pig, Cow, Cat, Dog |
| 免疫原 | A synthetic peptide corresponding to a region within the C terminal amino acids 1008-1057 (LPLLATDPSS FYSIEKDGYN NNIHCLTKAI IQVSAALFTL YNKNIETHLK) of human NCKAP1L (NP_005328). Run BLAST with ExPASy Run BLAST with NCBI |
| 阳性对照 | 293T cell lysate |

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. |
| 存储溶液 | Preservative: None Constituents: 2% Sucrose, PBS |
| 纯度 | Immunogen affinity purified |
| 纯化说明 | Purified by peptide affinity chromatography method. |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab82664** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|-------|------|---|
| WB | | Use a concentration of 1 µg/ml. Predicted molecular weight: 128 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T. |
| ELISA | | Use at an assay dependent concentration. ELISA titre using peptide based assay, 1:62500. |

靶标

| | |
|-------|--|
| 组织特异性 | Expressed only in cells of hematopoietic origin. |
| 序列相似性 | Belongs to the HEM-1/HEM-2 family. |
| 细胞定位 | Cell membrane. |

图片



Anti-NCKAP1L antibody (ab82664) at 1 µg/ml
(in 5% skim milk / PBS buffer) + 293T cell
lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG diluted in
1/50000 - 1/100000 as secondary antibody.

Predicted band size: 128 kDa

Observed band size: 128 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors