

Anti-NGDN antibody ab124128

1 图像

概述

| | |
|-------|--|
| 产品名称 | Anti-NGDN抗体 |
| 描述 | 兔多克隆抗体to NGDN |
| 经测试应用 | 适用于: WB |
| 种属反应性 | 与反应: Human |
| 免疫原 | KLH conjugated synthetic peptide designed between amino acids 135-165 from the Central region of Human NGDN (NP_056329.1, NP_001036100.1). |
| 阳性对照 | HepG2 cell lysate |

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term. |
| 存储溶液 | Preservative: 0.09% Sodium azide Constituent: 99% PBS |
| 纯度 | Immunogen affinity purified |
| 纯化说明 | This antibody is purified through a protein A column, followed by peptide affinity purification. |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab124128** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----|------|--|
| WB | | 1/100 - 1/500. Predicted molecular weight: 36 kDa. |

靶标

| | |
|----|---|
| 功能 | Involved in the translational repression of cytoplasmic polyadenylation element (CPE)-containing mRNAs. |
|----|---|


序列相似性

Belongs to the SAS10 family.

细胞定位

Nucleus. Cytoplasm. Cell projection > axon. Cell projection > dendrite. Cell projection > filopodium. Detected in axons, dendrites and filopodia. Colocalized with EIF4E in neurites.

Anti-NGDN antibody 图像



Anti-NGDN antibody (ab124128) at 1/100 dilution + HepG2 cell line lysates at 35 µg

Predicted band size : 36 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors