

## Product datasheet

# Anti-NNT1 antibody (Biotin) ab83408

### 概述

产品名称	Anti-NNT1抗体(Biotin)
描述	兔多克隆抗体to NNT1 (Biotin)
宿主	Rabbit
偶联物	Biotin
经测试应用	适用于: WB, Sandwich ELISA
种属反应性	与反应: Human
免疫原	Highly pure (>98%) recombinant human NNT1

### 性能

形式	Lyophilised:Centrifuge vial prior to opening. Reconstitute in sterile PBS containing 0.1% BSA to a concentration of 0.1-1.0 mg/ml.
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: None Constituents: PBS, pH 7.2
纯化说明	ab83408 is affinity purified and sterile filtered.
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab83408** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 0.1 - 0.2 µg/ml. Predicted molecular weight: 25 kDa.
Sandwich ELISA		Use a concentration of 0.25 - 1 µg/ml.

## 靶标

---

<b>功能</b>	Cytokine with B-cell stimulating capability. Binds to and activates the ILST/gp130 receptor.
<b>组织特异性</b>	Expressed predominantly in lymph nodes, spleen, peripheral blood lymphocytes, bone marrow, and fetal liver.
<b>疾病相关</b>	Defects in CLCF1 are the cause of cold-induced sweating syndrome type 2 (CISS2) [MIM:610313]. Cold-induced sweating syndrome (CISS) is an autosomal recessive disorder characterized by profuse sweating induced by cool surroundings (temperatures of 7 to 18 degrees Celsius). Additional abnormalities include a high-arched palate, nasal voice, depressed nasal bridge, inability to fully extend the elbows and kyphoscoliosis.
<b>序列相似性</b>	Belongs to the IL-6 superfamily.
<b>细胞定位</b>	Secreted.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors