

Product datasheet

Natural Adenovirus type 40 protein ab123997

概述

产品名称	Natural Adenovirus type 40 protein
蛋白长度	Full length protein

描述

性质	Native
来源	Native

氨基酸序列

技术指标

Our [Abpromise guarantee](#) covers the use of **ab123997** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

形式	Liquid
----	--------

补充说明	Inactivation: Heated to 54 degrees for 2 hours. Confirmed inactive by growth on 293 cells for 2 weeks with no CPE observed.
------	-----------------------------------------------------------------------------------------------------------------------------

This protein is derived from a pathogenic organism, and may be involved in a disease process, consequently exposure may have adverse health effects.

制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. Constituents: 88% PBS, 1.19% HEPES, 10% Glycerol
--------	-----------------------------------------------------------------------------------------------------------------------------------------------

常规信息

相关性	Adenoviruses are DNA viruses that cause upper respiratory tract infections, conjunctivitis, and other infections in humans. Adenoviruses were discovered in 1953. About 47 different types have been identified since then, and about half of them are believed to cause human diseases. Infants and children are most commonly affected by adenoviruses. Adenovirus infections can occur throughout the year, but seem to be most common from fall to spring. Adenoviruses are responsible for 3-5% of acute respiratory infections in children and 2% of respiratory illnesses in
-----	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

adults.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors