

Product datasheet

Anti-Neurotrophin 3 antibody ab67484

概述

产品名称	Anti-Neurotrophin 3抗体
描述	羊多克隆抗体to Neurotrophin 3
宿主	Sheep
特异性	Less than 1% cross reactivity against mouse NGF and recombinant human BDNF. Approximately 5% cross-reactivity against Neurotrophin 4.
经测试应用	适用于: IHC-G, ELISA, Inhibition Assay, WB
种属反应性	与反应: Rat, Chicken, Human
免疫原	Recombinant full length protein (Human)

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: PBS
纯度	Protein G purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab67484** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-G		Use at an assay dependent dilution.
ELISA		Use a concentration of 1 - 10 µg/ml.
Inhibition Assay		Use a concentration of 1 - 10 µg/ml.

应用	Ab评论	说明
WB		Use a concentration of 1 - 10 µg/ml. Predicted molecular weight: 29 kDa.

靶标

功能	Seems to promote the survival of visceral and proprioceptive sensory neurons.
组织特异性	Brain and peripheral tissues.
序列相似性	Belongs to the NGF-beta family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors