

## Product datasheet

# Anti-OSCP antibody ab93592

### 1 图像

#### 概述

产品名称	Anti-OSCP抗体
描述	兔多克隆抗体to OSCP
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human
免疫原	Synthetic peptide selected from the central region of Human OSCP conjugated to KLH (NP_004782.1).
阳性对照	Hela cell line lysates.

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Immunogen affinity purified
纯化说明	ab93592 is purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab93592** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB

ELISA

应用说明	ELISA: 1/1000. WB: 1/100 - 1/500. Predicted molecular weight: 54 kDa.
------	--


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## 靶标

相关性	Osteoblast specific cysteine-rich protein (OSCP) is integrin, $\beta$ -like 1 (with EGF-like repeat domains).
细胞定位	Secreted

## 图片



Anti-OSCP antibody (ab93592) at 1/100 dilution + Hela cell line lysates at 35  $\mu$ g

**Predicted band size : 54 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors