


Product datasheet

Anti-PAQR4 antibody ab125435

1 图像

概述

产品名称	Anti-PAQR4抗体
描述	兔多克隆抗体to PAQR4
经测试应用	适用于: WB
种属反应性	与反应: Rat 预测可用于: Mouse, Rabbit, Goat, Guinea pig, Cow, Dog, Human 
免疫原	Synthetic peptide corresponding to a region within C terminal amino acids 140-189 (ACRPWLRPAA LMGYTALSGV AGWRALTAPS TSARLRAFGW QAGARLLVFG) of Rat PAQR4 (NP_001017377; UniProtKB/TrEMBL ID: Q568Z3). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Rat lung lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituents: 97% PBS, 2% Sucrose
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用


Our [Abpromise guarantee](#) covers the use of **ab125435** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 29 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

组织特异性	Relatively widely expressed in a range of tissues.
序列相似性	Belongs to the ADIPOR family.
细胞定位	Membrane.

Anti-PAQR4 antibody 图像



Anti-PAQR4 antibody (ab125435) at 1 μ g/ml +
Rat lung lysate at 10 μ g

Predicted band size : 29 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors