

## Product datasheet

# Anti-PARVB antibody ab80691

★★★★☆ 1 Abreviews 1 图像

### 概述

产品名称	Anti-PARVB抗体
描述	兔多克隆抗体to PARVB
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Zebrafish
免疫原	Synthetic peptide corresponding to Human PARVB aa 41-90 (N terminal). Sequence: LQEEGKNAINSPMSPALVDVHPEDTQLEENEERTMIDPTSKEDPKFKELV <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
阳性对照	Transfected 293T cell lysate.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab80691** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
----	------	----

WB ★★★★★ Use a concentration of 1 µg/ml. Detects a band of approximately 45 kDa (predicted molecular weight: 42 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

功能	Probably plays a role in the regulation of cell adhesion and cytoskeleton organization.
组织特异性	Expressed predominantly in heart and skeletal muscle.
序列相似性	Belongs to the parvin family. Contains 2 CH (calponin-homology) domains.
细胞定位	Cell junction > focal adhesion. Cell membrane. Cytoplasm > cytoskeleton. Constituent of focal adhesions.

## 图片



Anti-PARVB antibody (ab80691) at 1 µg/ml + transfected 293T cell lysate at 10 µg

### Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

**Predicted band size:** 42 kDa

**Observed band size:** 45 kDa

Gel concentration 12%

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors