


Product datasheet

Anti-PDIR antibody ab97328

1 图像

概述

产品名称	Anti-PDIR抗体
描述	兔多克隆抗体to PDIR
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Cow 
免疫原	Recombinant fragment corresponding to a region within amino acids 1-223 of Human PDIR (AAH01625; UniProt ID: Q14554).
阳性对照	HepG2 whole cell lysate; A549 and HCT116 cell lysates

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 20% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab97328** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 60 kDa.

## 靶标


### 序列相似性

Belongs to the protein disulfide isomerase family.  
Contains 3 thioredoxin domains.

### 细胞定位

Endoplasmic reticulum lumen.

## 图片



Anti-PDIR antibody (ab97328) at 1/1000  
dilution + HepG2 whole cell lysate at 30  $\mu$ g

**Predicted band size:** 60 kDa

12% SDS-PAGE

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors