

Product datasheet

Anti-PISD antibody ab93603

1 图像

概述

产品名称	Anti-PISD抗体
描述	兔多克隆抗体to PISD
宿主	Rabbit
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Mouse, Human
免疫原	KLH conjugated synthetic peptide selected from the center region of Human PISD (NP_055153.1)
阳性对照	NIH3T3 cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Immunogen affinity purified
纯化说明	This antibody is purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab93603** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		
ELISA		

应用说明

ELISA: 1/1000.

WB: 1/100 - 1/500. Predicted molecular weight: 47 kDa.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

通路

Phospholipid metabolism; phosphatidylethanolamine biosynthesis; phosphatidylethanolamine from CDP-diacylglycerol: step 2/2.


序列相似性

Belongs to the phosphatidylserine decarboxylase family.

细胞定位

Mitochondrion.

图片



Anti-PISD antibody (ab93603) at 1/100 dilution + NIH3T3 cell lysate at 35 µg

Predicted band size: 47 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors