


Product datasheet

Anti-PLA2G1B antibody ab125866

1 图像

概述

产品名称	Anti-PLA2G1B抗体
描述	兔多克隆抗体to PLA2G1B
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Dog 
免疫原	Recombinant fragment, corresponding to a region within N terminal amino acids 1-136 of Human PLA2G1B (NP_000919).
阳性对照	A549 and HeLa whole cell lysates.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 0.75% Glycine, 1.21% Tris, 20% Glycerol, 1% BSA
纯度	Immunogen affinity purified
纯化说明	ab125866 is purified by antigen-affinity chromatography.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab125866** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 16 kDa.

靶标

功能	PA2 catalyzes the calcium-dependent hydrolysis of the 2-acyl groups in 3-sn-phosphoglycerides.
序列相似性	Belongs to the phospholipase A2 family.
细胞定位	Secreted.

图片



Anti-PLA2G1B antibody (ab125866) at 1/500 dilution + HeLa whole cell lysate at 30 µg

Predicted band size: 16 kDa

15% SDS PAGE

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors