

## Product datasheet

# Anti-PP4R1 antibody ab70624

### 2 图像

#### 概述

产品名称	Anti-PP4R1抗体
描述	兔多克隆抗体to PP4R1
宿主	Rabbit
经测试应用	适用于: WB, IP
种属反应性	与反应: Mouse, Human 预测可用于: Rhesus monkey
免疫原	Synthetic peptide, corresponding to amino acids 1-50 of Human PP4R1 (NP_005125.1)
阳性对照	HeLa whole cell lysate 293T whole cell lysate NIH3T3 whole cell lysate

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: 0.1% BSA, Tris buffered saline
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab70624** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/2000 - 1/10000. Detects a band of approximately 117 kDa (predicted molecular weight: 107 kDa).
IP		Use at 2-5 µg/mg of lysate.

## 靶标


## 相关性

PP4R1 is a protein which is associated with the catalytic subunit of Protein Serine/Threonine Phosphatase 4 (PP4C) and may possibly regulate the activity of the enzyme.

## 细胞定位

Nuclear and cytoplasmic.

## 图片



**All lanes :** Anti-PP4R1 antibody (ab70624) at 0.04 µg/ml

**Lane 1 :** HeLa whole cell lysate at 50 µg

**Lane 2 :** HeLa whole cell lysate at 15 µg

**Lane 3 :** HeLa whole cell lysate at 5 µg

**Lane 4 :** 293T whole cell lysate at 50 µg


**Lane 5 :** NIH3T3 whole cell lysate at 50 µg

**Predicted band size:** 107 kDa

**Observed band size:** 117 kDa

**Additional bands at:** 100 kDa, 120 kDa, 130 kDa, 70 kDa, 95 kDa. We are unsure as to the identity of these extra bands.

**Exposure time:** 3 minutes



HeLa whole cell lysate (1mg) was immunoprecipitated using ab70624 at 3µg/mg (lane 1) or control IgG (lane 2). 20% of the immunoprecipitate was used for Western blot and the bands revealed using ab70624 at 1µg/ml. Chemiluminescence was used for detection with an exposure time of 10 seconds.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors