

Product datasheet

Anti-PR3 antibody ab88097

2 图像

概述

产品名称	Anti-PR3抗体
描述	小鼠多克隆抗体to PR3
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length human PR3 protein (AAH96184.1).
阳性对照	Human spleen tissue lysate and PR3 transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab88097** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 28 kDa.

靶标

<b>功能</b>	Polymorphonuclear leukocyte serine protease that degrades elastin, fibronectin, laminin, vitronectin, and collagen types I, III, and IV (in vitro) and causes emphysema when administered by tracheal insufflation to hamsters.
<b>序列相似性</b>	Belongs to the peptidase S1 family. Elastase subfamily. Contains 1 peptidase S1 domain.


## 图片



Anti-PR3 antibody (ab88097) at 1  $\mu\text{g/ml}$  +  
hman spleen tissue lysate at 50  $\mu\text{g}$

**Predicted band size:** 28 kDa

**Observed band size:** 25 kDa



**All lanes :** Anti-PR3 antibody (ab88097) at 1  $\mu\text{g/ml}$

**Lane 1 :** PR3 transfected 293T cell lysate

**Lane 2 :** Non transfected 293T cell lysate

Lysates/proteins at 50  $\mu\text{g}$  per lane.

**Predicted band size:** 28 kDa

**Observed band size:** 31 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors