

Product datasheet

Anti-PTOV1 antibody ab106964

1 图像

概述

产品名称	Anti-PTOV1抗体
描述	兔多克隆抗体to PTOV1
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse, Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to N terminal amino acids 28-56 of Human PTOV1 (NP_059128.2).
阳性对照	mouse cerebellum tissue lysates

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Immunogen affinity purified
纯化说明	This antibody is purified through a protein A column, followed by peptide affinity purification.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab106964** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/100 - 1/500. Predicted molecular weight: 47 kDa.

## 靶标


### 相关性

PTOV1 belongs to the Mediator complex subunit 25 family, PTOV1 subfamily. It may activate transcription and is required for nuclear translocation of FLOT1. PTOV1 promotes cell proliferation. There are three named isoforms.

### 细胞定位

Cell membrane. Cytoplasm, perinuclear region. Nucleus. Translocates from the cytoplasm to the nucleus at the onset of S-phase. Also localizes to lipid rafts.

## 图片



Anti-PTOV1 antibody (ab106964) at 1/100 dilution + mouse cerebellum tissue lysate at 35  $\mu$ g

**Predicted band size:** 47 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors