

Product datasheet

Anti-Paramyxovirus Type 3 antibody ab33637

概述

产品名称	Anti-Paramyxovirus Type 3抗体
描述	鸡多克隆抗体to Paramyxovirus Type 3
宿主	Chicken
特异性	This antibody is monospecific.
经测试应用	适用于: ELISA, Neutralising
种属反应性	Reacts with Paramyxovirus.
免疫原	Paramyxovirus Type 3.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Add glycerol to a final volume of 50% for extra stability and aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: None Constituents: Whole serum
纯度	Whole antiserum
纯化说明	This antibody is a sterile filtered serum produced in SPF (Specific Pathogen Free) chicken embryos.
克隆	多克隆
同种型	IgY

应用

Our [Abpromise guarantee](#) covers the use of **ab33637** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.
Neutralising		Use at an assay dependent dilution.

靶标

相关性

The family of Paramyxoviridae contains viruses that induce a wide range of distinct clinical illnesses in humans. These include measles virus, which in rare instances is followed by subacute sclerosing panencephalitis (SSPE); mumps virus, which has symptoms of parotitis, orchitis and encephalitis, and the parainfluenza viruses which are respiratory pathogens.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors