

## Product datasheet

# Anti-Prostate Specific Antigen antibody ab19554

### 1 References

#### 概述

产品名称	Anti-Prostate Specific Antigen抗体
描述	兔多克隆抗体to Prostate Specific Antigen
经测试应用	适用于: ELISA
种属反应性	与反应: Human
免疫原	Full length protein (Human).

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: 0.85% Sodium chloride, 0.015M Potassium phosphate. pH 7.4
纯度	IgG fraction
纯化说明	DEAE column chromatography.
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab19554** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

ELISA

应用说明	ELISA: Suitable for solid phase coating. Use at an assay dependent dilution.  Not tested in other applications. Optimal dilutions/concentrations should be determined by the end user.
------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

## 靶标

---

功能	Hydrolyzes semenogelin-1 thus leading to the liquefaction of the seminal coagulum.
序列相似性	Belongs to the peptidase S1 family. Kallikrein subfamily. Contains 1 peptidase S1 domain.
细胞定位	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors