

Product datasheet

Anti-RASA3 antibody ab31286


概述

产品名称

Anti-RASA3抗体

描述

山羊多克隆抗体to RASA3

 该产品是 **快线抗体**，通过亲和纯化得到。对抗原肽有高灵敏度及特异性 (ELISA测试)。

[Read the terms of use »](#)

宿主

Goat

种属反应性

预测可用于: Mouse, Rat, Dog, Human 

免疫原

Synthetic peptide:

C-QYKRDKFKKTKYGSQ

, corresponding to C terminal amino acids 796-810 of Human RASA3

 [Run BLAST with](#)

 [Run BLAST with](#)

性能

形式

Liquid

存放说明

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

存储溶液

Preservative: 0.02% Sodium Azide

Constituents: 0.5% BSA, Tris buffered saline, pH 7.3

纯度

Immunogen affinity purified

纯化说明

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

克隆

多克隆

同种型

IgG

应用

应用说明

This antibody gave a positive result in ELISA against the immunizing peptide. Antibody detection limit dilution 1:32,000.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

---

## 靶标

功能	Inhibitory regulator of the Ras-cyclic AMP pathway. Binds inositol tetrakisphosphate (IP4) with high affinity. Might be a specific IP4 receptor.
序列相似性	Contains 1 Btk-type zinc finger. Contains 2 C2 domains. Contains 1 PH domain. Contains 1 Ras-GAP domain.
细胞定位	Cell membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

---

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors