

## Product datasheet

# Anti-RCN2 antibody ab104516

### 2 图像

#### 概述

产品名称	Anti-RCN2抗体
描述	兔多克隆抗体to RCN2
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length Human RCN2 (NP_002893.1; amino acids 1-317).
阳性对照	HeLa cell lysate; RCN2 transfected 293T cell lysate

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab104516** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 37 kDa.

#### 靶标

功能	Not known. Binds calcium.
----	---------------------------


组织特异性	Ubiquitous.
序列相似性	Belongs to the CREC family. Contains 6 EF-hand domains.
细胞定位	Endoplasmic reticulum lumen.

## 图片



Anti-RCN2 antibody (ab104516) at 1  $\mu$ g/ml +  
HeLa cell lysate at 50  $\mu$ g

**Predicted band size:** 37 kDa



**All lanes :** Anti-RCN2 antibody (ab104516) at  
1  $\mu$ g/ml

**Lane 1 :** RCN2 transfected 293T cell lysate

**Lane 2 :** Non-transfected 293T cell lysate

Lysates/proteins at 25  $\mu$ g per lane.

**Predicted band size:** 37 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors