

Product datasheet

Rabbit Anti-Goat IgG Fc (DyLight® 650) preadsorbed ab102343

概述

产品名称	兔抗山羊IgG Fc (DyLight® 650)预吸附二抗
宿主	Rabbit
靶标种属	Goat
特异性	By immunoelectrophoresis and ELISA this antibody reacts specifically with goat IgG. No antibody was detected against non-immunoglobulin serum proteins. Less than 1% cross reactivity to human, mouse, rat IgG was detected. This antibody may cross react with IgG from other species.
经测试应用	适用于: IHC-P, ICC, Flow Cyt
最小交叉反应	Human, Mouse, Rat more details
偶联物	DyLight® 650. Ex: 654nm, Em: 673nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Immunogen affinity purified
纯化说明	Antiserum was solid phase absorbed to ensure class specificity. Antiserum was cross absorbed using human, mouse, rat immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. Fragments were conjugated to DyLight® 650.
克隆	多克隆
同种型	IgG
常规说明	DyLight® 650 replaces DyLight® 649.

应用

Our [Abpromise guarantee](#) covers the use of **ab102343** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		1/50 - 1/500.
ICC		1/50 - 1/500.
Flow Cyt		1/50 - 1/200.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors