

## Product datasheet

# Rabbit Anti-Goat IgG Fc (DyLight® 650) preadsorbed ab102343

### 概述

<b>产品名称</b>	兔Anti-Goat IgG Fc (DyLight® 650) preadsorbed
<b>靶标种属</b>	Goat
<b>特异性</b>	By immunoelectrophoresis and ELISA this antibody reacts specifically with goat IgG. No antibody was detected against non-immunoglobulin serum proteins. Less than 1% cross reactivity to human, mouse, rat IgG was detected. This antibody may cross react with IgG from other species.
<b>经测试应用</b>	<b>适用于:</b> IHC-P, ICC, Flow Cyt
<b>最小交叉反应</b>	Human, Mouse, Rat <a href="#">more details</a>
<b>偶联物</b>	DyLight® 650. Ex: 654nm, Em: 673nm

### 性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C.
<b>存储溶液</b>	Preservative: 0.09% Sodium Azide Constituents: PBS
<b>纯度</b>	Immunogen affinity purified
<b>纯化说明</b>	Antiserum was solid phase absorbed to ensure class specificity. Antiserum was cross absorbed using human, mouse, rat immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. Fragments were conjugated to DyLight® 650.
<b>克隆</b>	多克隆
<b>同种型</b>	IgG
<b>常规说明</b>	DyLight® 650 replaces DyLight® 649.

### 应用

Our [Abpromise guarantee](#) covers the use of **ab102343** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		1/50 - 1/500.
ICC		1/50 - 1/500.
Flow Cyt		1/50 - 1/200.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors