

## Product datasheet

# Rabbit Anti-Guinea pig IgG H&L (FITC) ab6767

## 2 References

### 概述

产品名称	兔Anti-Guinea猪IgG H&L (FITC)
靶标种属	Guinea pig
经测试应用	适用于: IHC-P, Immunomicroscopy, Flow Cyt, ICC/IF, IHC-Fr, ELISA
免疫原	guinea pig whole molecule
偶联物	FITC. Ex: 493nm, Em: 528nm

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.2 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 0.164% Sodium phosphate, 0.878% Sodium chloride, 1% BSA
纯度	Affinity purified
纯化说明	Diaminotriazinylaminofluorescein (DTAF) (Molecular Weight 530 daltons) Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Fluorochrome/Protein Ratio: 2.3 moles DTAF per mole of Rabbit IgG
共轭说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using guinea pig IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
克隆	多克隆

### 应用

Our [Abpromise guarantee](#) covers the use of **ab6767** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		Use at an assay dependent concentration.
Immunomicroscopy		Use at an assay dependent concentration.

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent dilution.
ICC/IF		1/1000 - 1/5000.
IHC-Fr		Use at an assay dependent concentration.
ELISA		1/10000 - 1/50000.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors