

Product datasheet

Rabbit Anti-Guinea pig IgG H&L (FITC) ab6767

2 References

概述

| | |
|-------|---|
| 产品名称 | 兔Anti-Guinea猪IgG H&L (FITC) |
| 靶标种属 | Guinea pig |
| 经测试应用 | 适用于: IHC-P, Immunomicroscopy, Flow Cyt, ICC/IF, IHC-Fr, ELISA |
| 免疫原 | guinea pig whole molecule |
| 偶联物 | FITC. Ex: 493nm, Em: 528nm |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C. |
| 存储溶液 | pH: 7.2 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 0.164% Sodium phosphate, 0.878% Sodium chloride, 1% BSA |
| 纯度 | IgG fraction |
| 纯化说明 | Diaminotriazinylaminofluorescein (DTAF) (Molecular Weight 530 daltons) Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Fluorochrome/Protein Ratio: 2.3 moles DTAF per mole of Rabbit IgG |
| 共轭说明 | This product was prepared from monospecific antiserum by immunoaffinity chromatography using guinea pig IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. |
| 克隆 | 多克隆 |

应用

Our [Abpromise guarantee](#) covers the use of **ab6767** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|------------------|------|--|
| IHC-P | | Use at an assay dependent concentration. |
| Immunomicroscopy | | Use at an assay dependent concentration. |

| 应用 | Ab评论 | 说明 |
|----------|------|--|
| Flow Cyt | | Use at an assay dependent dilution. |
| ICC/IF | | 1/1000 - 1/5000. |
| IHC-Fr | | Use at an assay dependent concentration. |
| ELISA | | 1/10000 - 1/50000. |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors