

## Product datasheet

# Rabbit Anti-Horse IgG+IgM+IgA H&L (HRP) ab8533

### 概述

产品名称	兔抗马IgG+IgM+IgA H&L (HRP)
宿主	Rabbit
靶标种属	Horse
特异性	Antibodies were raised in rabbit to Horse IgG, IgM and IgA
经测试应用	适用于: ELISA, Dot blot
偶联物	HRP

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	50% glycerol, 0.1 M NaCl, 0.01 M sodium phosphate at pH 7.5.
纯度	Immunogen affinity purified
纯化说明	Specific antibodies were absorbed by incubation with Sepharose-bound Horse immunoglobulins. Specific antibodies were eluted by acidic buffer at pH 2.5 followed by neutralisation and dialysis. After repeated binding with immobilized Horse immunoglobulins a minimum of 65% protein bound.
共轭说明	The purified polyclonal was conjugated to horseradish peroxidase according to the periodate method followed by gel-filtration and ion-exchange chromatography to clear unbound conjugate.
克隆	多克隆

### 应用

Our [Abpromise guarantee](#) covers the use of **ab8533** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

ELISA

Dot blot

**应用说明** Recommended dilution as secondary antibody in ELISA is 1:1000.

Recommended dilution as secondary antibody in dot-blotting is 1:100.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors