

Product datasheet

Rabbit Anti-Human IgM mu chain ab97206

概述

| | |
|-------|--|
| 产品名称 | 兔抗人IgM mu chain |
| 宿主 | Rabbit |
| 靶标种属 | Human |
| 特异性 | By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgM. Cross reactivity with other immunoglobulins and light chains is less than 0.1%. |
| 经测试应用 | 适用于: ICC, IHC-P, ELISA, WB |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| 存储溶液 | Preservative: 0.09% Sodium Azide Constituents: PBS |
| 纯度 | Immunogen affinity purified |
| 纯化说明 | Antiserum was solid phase adsorbed to ensure class specificity. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab97206** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|-------|------|--|
| ICC | | Use at an assay dependent dilution. |
| IHC-P | | 1/200 - 1/2000. |
| ELISA | | 1/100 - 1/30000. ELISA-Coating 1/100 - 1/500; ELISA-Primary 1/1000 - 1/30000 |

| 应用 | Ab评论 | 说明 |
|----|------|--|
| WB | | 1/1000 - 1/30000. WB-Colorimetric 1/1000 - 1/10000; WB-Chemiluminescent 1/1000 - 1/30000 |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors