

## Product datasheet

# Rabbit Anti-Pig IgG H&L (TRITC) ab6774

### 概述

产品名称	兔Anti-Pig IgG H&L (TRITC)
靶标种属	Pig
经测试应用	<b>适用于:</b> Immunomicroscopy, Flow Cyt, ICC/IF, IHC-P, IHC-Fr, ELISA
免疫原	pig IgG whole molecule
偶联物	TRITC. Ex: 547nm, Em: 572nm

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.2 Preservative: 0.01% Sodium azide Constituents: 0.878% Sodium chloride, 1% BSA, 0.424% Potassium phosphate
纯度	IgG fraction
纯化说明	Tetramethylrhodamine isothiocyanate (TRITC) (Molecular Weight 444 daltons) Absorption Wavelength: 550 nm Emission Wavelength: 570 nm Fluorochrome/Protein Ratio: 2.8 moles TRITC per mole of Rabbit IgG
共轭说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using pig IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities
克隆	多克隆

### 应用

Our [Abpromise guarantee](#) covers the use of **ab6774** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Immunomicroscopy		
Flow Cyt		
ICC/IF		

---

IHC-P

---

IHC-Fr

---

ELISA

---

**应用说明**

Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring lot-to-lot consistency

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

---

**Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

---

**Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors