

## Product datasheet

# Rabbit Anti-Rat IgG H&L preadsorbed ab102248

### 1 References

#### 概述

<b>产品名称</b>	兔Anti-Rat IgG H&L preadsorbed
<b>描述</b>	兔多克隆抗体二抗to大鼠IgG - H&L, pre-adsorbed
<b>靶标种属</b>	Rat
<b>特异性</b>	By immunoelectrophoresis and ELISA this antibody reacts specifically with rat IgG and with light chains common to other rat immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. Less than 1% cross reactivity to bovine, horse, human, mouse, rabbit and sheep IgG was detected. This antibody may cross react with IgG from other species.
<b>经测试应用</b>	<b>适用于:</b> ICC, IHC-P, ELISA, WB
<b>最小交叉反应</b>	Cow, Horse, Human, Mouse, Rabbit, Sheep <a href="#">more details</a>

#### 性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>存储溶液</b>	Preservative: 0.09% Sodium Azide Constituents: PBS
<b>纯度</b>	Immunogen affinity purified
<b>纯化说明</b>	Antiserum was cross adsorbed using bovine, horse, human, mouse, rat and sheep immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads.
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab102248** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC		1/200 - 1/2000.

IHC-P	1/200 - 1/2000.
ELISA	1/100 - 1/30000. Coating: Use 1/100 - 1/500 dilution; Primary: Use 1/1000 - 1/30000 dilution.
WB	1/1000 - 1/30000. Colorimetric: Use 1/1000 - 1/10000 dilution; Chemiluminescent: Use 1/1000 - 1/30000 dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors