

Product datasheet

Rabbit Anti-Rat IgG H&L preadsorbed ab102248

1 References

概述

产品名称	兔抗大鼠IgG H&L预吸附二抗
描述	兔多克隆抗体二抗to大鼠IgG - H&L,预吸附二抗
宿主	Rabbit
靶标种属	Rat
特异性	By immunoelectrophoresis and ELISA this antibody reacts specifically with rat IgG and with light chains common to other rat immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. Less than 1% cross reactivity to bovine, horse, human, mouse, rabbit and sheep IgG was detected. This antibody may cross react with IgG from other species.
经测试应用	适用于: ICC, IHC-P, ELISA, WB
最小交叉反应	Cow, Horse, Human, Mouse, Rabbit, Sheep more details

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Immunogen affinity purified
纯化说明	Antiserum was cross adsorbed using bovine, horse, human, mouse, rat and sheep immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab102248** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC		1/200 - 1/2000.
IHC-P		1/200 - 1/2000.
ELISA		1/100 - 1/30000. Coating: Use 1/100 - 1/500 dilution; Primary: Use 1/1000 - 1/30000 dilution.
WB		1/1000 - 1/30000. Colorimetric: Use 1/1000 - 1/10000 dilution; Chemiluminescent: Use 1/1000 - 1/30000 dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors