

Product datasheet

Anti-Raspberry Bushy Dwarf Virus Coat Protein antibody ab12243

概述

产品名称	Anti-Raspberry Bushy Dwarf Virus Coat蛋白抗体
描述	鸡多克隆抗体to Raspberry Bushy Dwarf Virus Coat蛋白
宿主	Chicken
经测试应用	适用于: WB, ELISA
种属反应性	This antibody reacts with Raspberry Bushy Dwarf Virus coat protein of 32 kDa.
免疫原	Native Raspberry Bushy Dwarf Virus (RBDV) particles purified from Chenopodium quinoa.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.05% Sodium azide Constituents: PBS, 0.0292% EDTA
纯度	IgY fraction
克隆	多克隆
同种型	IgY

应用

Our [Abpromise guarantee](#) covers the use of **ab12243** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB

ELISA

应用说明	ELISA: Use at an assay dependent dilution. ELISA titer > 1:70,000. WB: Use at an assay dependent dilution. Detects a band of approximately 32 kDa.
------	---

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

相关性

Raspberry bushy dwarf virus, which is spread by pollen and infects many red and black raspberry cultivars. It is also found in 'Boysenberry', 'Loganberry' and 'Marionberry' blackberries. The fruit from infected plants often will not make IQF grade and is sold as a lower value produce for juice, jam or puree. Pollen is infected with the virus and honeybees move infected pollen from plant to plant.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors