

Product datasheet

Rat IgG (FITC) - Isotype Control ab37364

概述

| | |
|-------|---|
| 产品名称 | 大鼠IgG (FITC) -同型对照 |
| 偶联物 | FITC. Ex: 493nm, Em: 528nm |
| 特异性 | This product has been selected as an isotype control as it has no known specificity |
| 经测试应用 | 适用于: Flow Cyt, ELISA |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: 0.1% Sodium Azide Constituents: PBS |
| 纯度 | Affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab37364** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----------|------|-------------------------------------|
| Flow Cyt | | Use 1µg for 10 ⁶ cells. |
| ELISA | | Use at an assay dependent dilution. |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors