

Product datasheet

Rat P2X3 peptide ab46989

概述

产品名称 大鼠P2X3肽

描述

性质 Synthetic
 氨基酸序列
 种属 Rat
 序列 VEKQSTDSGAYSIGH
 氨基酸 383 to 397

技术指标

Our [Abpromise guarantee](#) covers the use of **ab46989** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 Blocking - Blocking peptide for Anti-P2X3 antibody ([ab10269](#))
 形式 Liquid

制备和贮存

稳定性和存储 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
 Peptide reconstituted in distilled water.

常规信息

功能 Receptor for ATP that acts as a ligand-gated ion channel.
 序列相似性 Belongs to the P2X receptor family.
 细胞定位 Membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors