

Product datasheet

Rat monoclonal [JC5-1] Anti-Mouse lambda light chain (Alkaline Phosphatase) ab99624

概述

产品名称	大鼠单克隆抗体[JC5-1]抗小鼠lambda light chain (Alkaline Phosphatase)
宿主	Rat
靶标种属	Mouse
特异性	ab99624 reacts with lambda 1 and lambda 2 light chains. It does not react with other immunoglobulin isotypes or kappa light chains.
经测试应用	适用于: ELISA, Flow Cyt, ICC/IF
偶联物	Alkaline Phosphatase

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: 50% Glycerol, 50mM Tris, 1mM Magnesium chloride, pH 8
纯度	Protein A purified
克隆	单克隆
克隆编号	JC5-1
同种型	IgG2b
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab99624** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		1/500 - 1/1000.
Flow Cyt		Use at an assay dependent dilution.

应用	Ab评论	说明
ICC/IF		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors