

## Product datasheet

# Rat monoclonal [LO-MG3-13] Anti-Mouse IgG3 heavy chain ab11656

### 概述

---

产品名称	大鼠单克隆抗体[LO-MG3-13]抗小鼠IgG3 heavy chain
宿主	Rat
靶标种属	Mouse
特异性	ab11656 recognises murine IgG3 (gamma heavy chain) and does not cross react with other murine immunoglobulin classes or subclasses.
经测试应用	适用于: ELISA
免疫原	Polyclonal mouse IgG.

### 性能

---

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Protein G purified
纯化说明	Purified IgG prepared by affinity chromatography.
克隆	单克隆
克隆编号	LO-MG3-13
骨髓瘤	IR983F
同种型	IgG1
常规说明	This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

### 应用

---

Our [Abpromise guarantee](#) covers the use of **ab11656** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use a concentration of 5 µg/ml.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors