

Product datasheet

Recombinant Human ANGPTL3 protein ab108555

概述

产品名称	重组人ANGPTL3蛋白
蛋白长度	Protein fragment

描述

性质	Recombinant
来源	HEK 293 cells
氨基酸序列	
Accession	Q9Y5C1
种属	Human
分子量	35 kDa including tags
氨基酸	224 to 460
标签	DDDDK tag N-Terminus

技术指标

Our [Abpromise guarantee](#) covers the use of **ab108555** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
内毒素水平	< 0.100 Eu/μg
纯度	> 90 % SDS-PAGE. ab108555 is 0.2 μm filtered.
形式	Liquid
补充说明	Working aliquots are stable for up to 3 months when stored at -20°C.

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Preservative: None Constituents: PBS
--------	--

常规信息

组织特异性	Expressed principally in liver. Weakly expressed in kidney.
序列相似性	Contains 1 fibrinogen C-terminal domain.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors