

Product datasheet

Recombinant Human ARL15 protein ab113118

1 图像

概述

产品名称	Recombinant人ARL15 protein
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli

氨基酸序列

Accession	Q9NXU5
种属	Human
序列	MGSSHHHHHSSGLVPRGSHMSDLRITEAFLYMDYLCFRALCCKGPPPAR PEYDLVCIGLTGSGKTSLLSKLCSESPDNVVSTTGFSIKAVPFQNAILNV KELGGADNIRKYWSRYQGSQGVIFVLDSASSEDDLEAARNELHSALQHP QLCTLPFLILANHQDKPAARSVQEIKKYFELEPLARGKRWILQPCSLDDM DALKDSFSQLINLLEEKDHEAVRM
分子量	25 kDa including tags
氨基酸	1 to 204
标签	His tag N-Terminus

技术指标

Our [Abpromise guarantee](#) covers the use of **ab113118** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
纯度	> 85 % SDS-PAGE. ab113118 is purified using conventional chromatography techniques.
形式	Liquid

稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00


Constituents: 0.32% Tris HCl, 0.02% DTT, 10% Glycerol, 0.58% Sodium chloride

常规信息

序列相似性

Belongs to the small GTPase superfamily. Arf family.

图片



15% SDS-PAGE analysis of 3 µg ab113118.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors