

Product datasheet

Recombinant Human CCL27 protein ab93953

1 图像

概述

产品名称	Recombinant人CCL27 protein
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli

氨基酸序列

种属	Human
序列	MFLLPSTAC CTQLYRKPLS DKLLRKVIQV ELQEADGDCH LQAFVLHLAQ RSICIHQPQP SLSQWFEHQE RKLHGTLPKL NFGMLRKMG

技术指标

Our [Abpromise guarantee](#) covers the use of **ab93953** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Mass Spectrometry SDS-PAGE
纯度	> 95 % SDS-PAGE. ab93953 was purified using conventional chromatography techniques.
形式	Liquid
补充说明	Endotoxin Level: < 1.0 EU per 1 µg of protein (determined by LAL method)


制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. Preservative: None Constituents: 10% Glycerol, 10mM Sodium citrate, pH 3.5
--------	---

常规信息

功能	Chemotactic factor that attracts skin-associated memory T-lymphocytes. May play a role in mediating homing of lymphocytes to cutaneous sites. Binds to CCR10.
组织特异性	Testis, thymus, placenta, ovary and skin.
序列相似性	Belongs to the intercrine beta (chemokine CC) family.
细胞定位	Secreted.

图片



15% SDS-PAGE analysis of ab93953 (3 μ g)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors