

## Product datasheet

# Recombinant Human CD252 protein ab69496

### 概述

---

产品名称	Recombinant人CD252 protein
蛋白长度	Full length protein

### 描述

---

性质	Recombinant
来源	Insect cells

### 氨基酸序列

种属	Human
----	-------

### 技术指标

---

Our [Abpromise guarantee](#) covers the use of **ab69496** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
内毒素水平	< 0.100 Eu/ $\mu$ g
纯度	> 95 % SDS-PAGE. Greater than 97% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per $\mu$ g (1EU/ $\mu$ g).
形式	Lyophilised

### 制备和贮存

---

稳定性和存储	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. lot specific
复溶	Lot specific

### 常规信息

---

功能	Cytokine that binds to TNFRSF4. Co-stimulates T-cell proliferation and cytokine production.
疾病相关	Genetic variations in TNFSF4 influence susceptibility to systemic lupus erythematosus (SLE)

[MIM:152700]. SLE is a chronic, inflammatory and often febrile multisystemic disorder of connective tissue. It affects principally the skin, joints, kidneys and serosal membranes. It is thought to represent a failure of the regulatory mechanisms of the autoimmune system. Note=The upstream region of TNFSF4 contains a single risk haplotype for SLE, which is correlated with increased expression of both cell-surface TNFSF4 and TNFSF4 transcripts. Increased levels of TNFSF4 are thought to augment T cell-APC interaction and the functional consequences of T cell activation, thereby destabilizing peripheral tolerance.

**序列相似性**

Belongs to the tumor necrosis factor family.

**细胞定位**

Membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors