

Product datasheet

Recombinant Human CHD4 protein ab114276

1 图像

概述

产品名称	Recombinant人CHD4 protein
蛋白长度	Protein fragment

描述

性质	Recombinant
来源	Wheat germ

氨基酸序列

Accession	<a href="#">Q14839</a>
种属	Human
序列	ETEPKGAADVKEKVEEKSAIDLTPIVVEDKEEKKEEEEKKEVMLQNGETPK DLNDEKQKKNIKQRFMFNIADGGFTELHSLWQNEERAATVTKKTYEIIWH
分子量	37 kDa including tags
氨基酸	1632 to 1730

技术指标

Our [Abpromise guarantee](#) covers the use of **ab114276** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA SDS-PAGE Western blot
形式	Liquid
补充说明	Protein concentration is above or equal to 0.05 mg/ml. This protein is best used within three months from the date of receipt.

制备和贮存


稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00
--------	---

Constituents: 0.79% Tris HCl, 0.3% Glutathione

## 常规信息

<b>功能</b>	Component of the histone deacetylase NuRD complex which participates in the remodeling of chromatin by deacetylating histones.
<b>序列相似性</b>	Belongs to the SNF2/RAD54 helicase family. Contains 2 chromo domains. Contains 1 helicase ATP-binding domain. Contains 1 helicase C-terminal domain. Contains 2 PHD-type zinc fingers.
<b>翻译后修饰</b>	Phosphorylated upon DNA damage, probably by ATM or ATR.
<b>细胞定位</b>	Nucleus. Cytoplasm > cytoskeleton > centrosome. Associates with centrosomes in interphase.

## 图片



12.5% SDS-PAGE analysis of CHD4 protein (ab114276). Stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors