

Product datasheet

Recombinant Human Calcipressin 1 protein ab63837

1 图像

概述

产品名称	重组人Calcipressin 1蛋白
蛋白长度	Protein fragment

描述

性质	Recombinant
来源	Escherichia coli

氨基酸序列

种属	Human
序列	MKLYFAQTLH IGSSHLAPPN PDKQFLISPP ASPPVGWKQV EDATPVINYD LLYAISKLG GEKYELHAAT DTPSVVVHV CESDQEKEEE EEMERMRRPK PKIIQTRRPE YTPIHLS
氨基酸	1 to 117

技术指标

Our [Abpromise guarantee](#) covers the use of **ab63837** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
纯度	> 90 % SDS-PAGE. ab63837 was expressed in E.coli and purified by using conventional chromatography techniques.
形式	Liquid


制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Preservative: None Constituents: 0.1M Sodium chloride, 20mM Tris 7.5, 1mM DTT
--------	---

## 常规信息

功能	Inhibits calcineurin-dependent transcriptional responses by binding to the catalytic domain of calcineurin A. Could play a role during central nervous system development.
组织特异性	Highly expressed heart, brain and skeletal muscle. Also expressed in all other tissues.
序列相似性	Belongs to the RCAN family.

## 图片



SDS Page analysis of ab63837 (3 ug on 15% gel).

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors