

Product datasheet

Recombinant Human Calcipressin 1 protein ab63837

1 图像

概述

产品名称	Recombinant人Calcipressin 1 protein
蛋白长度	Protein fragment

描述

性质	Recombinant
来源	Escherichia coli

氨基酸序列

种属	Human
序列	MKLYFAQTLH IGSSHLAPPN PDKQFLISPP ASPPVGWKQV EDATPVINYD LLYAISKLG GEKYELHAAT DTPSVVVHV CESDQEKEEE EEMERMRRPK PKIIQTRRPE YTPIHLS
氨基酸	1 to 117

技术指标

Our [Abpromise guarantee](#) covers the use of **ab63837** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
纯度	> 90 % SDS-PAGE. ab63837 was expressed in E.coli and purified by using conventional chromatography techniques.
形式	Liquid


制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Preservative: None Constituents: 0.1M Sodium chloride, 20mM Tris 7.5, 1mM DTT
--------	---

常规信息

功能	Inhibits calcineurin-dependent transcriptional responses by binding to the catalytic domain of calcineurin A. Could play a role during central nervous system development.
组织特异性	Highly expressed heart, brain and skeletal muscle. Also expressed in all other tissues.
序列相似性	Belongs to the RCAN family.

图片



SDS Page analysis of ab63837 (3 ug on 15% gel).

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors