

Product datasheet

Recombinant Human Cyclophilin F protein ab40517

1 图像

概述

产品名称	Recombinant人Cyclophilin F protein
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli
氨基酸序列	
种属	Human

技术指标

Our [Abpromise guarantee](#) covers the use of **ab40517** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Mass Spectrometry SDS-PAGE
纯度	> 95 % SDS-PAGE.
形式	Liquid

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. Preservative: 0.002% Sodium Azide Constituents: 10mM Tris, 10mM DTT, pH 8.0, 0.1% Triton X-100
--------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

常规信息

功能	PPlases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.
----	----------------------------------------------------------------------------------------------------------------------------------------

## 序列相似性


Belongs to the cyclophilin-type PPlase family.

Contains 1 PPlase cyclophilin-type domain.

## 细胞定位

Mitochondrion matrix.

## 图片



ab40517, 4-20% SDS gradient gel,  
coomassie blue staining. Expected size 22  
kDa.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors