

Product datasheet

Recombinant Human DNA polymerase delta p50 protein
ab114798

1 图像

概述

产品名称	重组人DNA polymerase delta p50蛋白
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Wheat germ

氨基酸序列

Accession [P49005](#)

种属 Human

序列 MFSEQAAQRAHTLLSPPSANNATFARVPVATYTNSSQPFR LGERSFSRQY
 AHIYATRLIQMRPFLENRAQQHWGSGVGVKKLCELQPEEKCCVGT LFKA
 MPLQPSILREVSEEHNLLQP PRSKYIHPDDELVLEDELQRIK LKGTIDV
 SKLVTGTVLAVFGSVRDDGKFLVEDYCFADLAPQKPAPPLDTRFVLLVS
 GLGLGGGGESLLGTQLLDV VVTGQLGDEGEQCSAAHVS RVILAGNLLSH
 STQSRDSINKAKYLTKKTQAASVEAVKMLDEILLQLSASVPVDVMPGEFD
 PTNYTLPQQPLHPCMFPLATAYSTLQLVTNPYQATIDGVRFLGTSGQNVS
 DIFRYSSMEDHLEILEWTLRVRHISPTAPDTLGCYPFYKTDPFIFPECPH
 VYFCGNTPSFGSKIIRGPEDQTVLLVTVPDFSATQTA CLVNLRSLACQPI
 SFSGFGAEDDDLGLGLGP

分子量 78 kDa including tags

氨基酸 1 to 469

技术指标

Our [Abpromise guarantee](#) covers the use of **ab114798** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot
	SDS-PAGE
	ELISA

形式	Liquid
补充说明	Protein concentration is above or equal to 0.05 mg/ml. Best used within three months from the date of receipt.


制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
--------	---------------------------------------------------------------------------------------------------------------------------------------------------------

常规信息

功能	The function of the small subunit is not yet clear.
序列相似性	Belongs to the DNA polymerase delta/II small subunit family.
细胞定位	Nucleus.

图片



ab114798 on 12.5% SDS-PAGE Stained with Coomassie Blue.

SDS-PAGE - DNA polymerase delta p50 protein
(ab114798)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors